

CM3215

MichiganTech

Fundamentals of Chemical Engineering Laboratory

Viscosity Measurement


Professor Faith Morrison

Department of Chemical Engineering
Michigan Technological University

1

© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory



We can use flow through a tube to measure viscosity:

Hagen-Poiseuille Equation:

without gravity $Q = \frac{\pi(p_0 - p_L)R^4}{8\mu L}$

with gravity $Q = \frac{\pi(p_0 - p_L + \rho g L \cos \beta)R^4}{8\mu L}$

Assumptions:

- Steady state
- Incompressible
- Newtonian
- No end effects

(when the capillary is tilted, is always constant)

Reference: Morrison, *An Introduction to Fluid Mechanics*, Cambridge, 2013.

2

© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory

Cannon-Fenske Viscometer

- Gravity driven
- Results in *kinematic viscosity*
- Requires loading a fixed volume that is specific to the supplied calibrations

$$v = \frac{\mu}{\rho}$$

3
© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory

$Q = \frac{\Delta V}{\Delta t}$ (volume)

$Q = \frac{\pi (p_0 - p_L + \rho g L \cos \beta) R^4}{8\mu L}$ Hagen-Poiseuille

Neglect Δp , solve for v :

$v = \frac{\mu}{\rho} = \left(\frac{\text{correction}}{\text{factor}} \right) \left(\frac{\pi R^4 g \cos \beta}{8\Delta V} \right) \Delta t$ (volume)

Kinematic Viscosity, v $v = \alpha \Delta t$

capillary tube

(For complete analysis, see Morrison, *An Introduction to Fluid Mechanics*, Example 7.4, p508)

$\mu = \rho \alpha \Delta t$

4
© Faith A. Morrison, Michigan Tech U.


CM3215 Fundamentals of Chemical Engineering Laboratory

Ubbelohde Viscometer

Kinematic Viscosity, ν

$$\nu = \frac{\mu}{\rho}$$

- Gravity driven
- Results in kinematic viscosity
- **Vented exit allows variable volumes to be charged**
- Design eliminates Δp




5

© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory

Loading a Cannon-Fenske Viscometer

- Accurate calibration requires using the calibrated volume
- Fill to top mark



suction

6


© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory

See class website for detailed instructions on operating the *Cannon-Fenske Routine Viscometers*

Instructions for the use of the Cannon-Fenske Routine Viscometer
 Cannon Instrument Company
 Modified by F. A. Morrison 13 September 2007

1. Obtain a stop watch, a viscometer of appropriate size, and approximately 30 ml of sample to be measured.
2. Clean the viscometer using suitable solvents and by passing clean, dry, filtered air through the instrument to remove the final traces of solvent (clean, dry air is available at two locations on the cast wall over the first-aid kit). Periodically, traces of organic deposits should be removed with chromic acid or non-chromium cleaning solution. Be careful when blowing air into the viscometer; do not drop.
3. If there is a possibility of lint, dust, or other solid material in the liquid sample, filter the sample through a sintered glass filter or fine mesh screen.
4. To charge the sample into the viscometer, invert the instrument and apply suction to tube G, immersing tube A in the liquid sample, and draw liquid to mark E (this can be accomplished by using the baby ear syring provided; these seal well on the inside of tube G). Turn the instrument to its normal vertical position and wipe clean arm A. This procedure charges the viscometer with the correct amount of liquid for the calibration constant given.
5. Place the viscometer into the constant temperature bath using a ring stand and a clamp. Align the tops of the tubes of the viscometer vertically in the bath. Up to four viscometers may be present in the bath at one time.
6. Allow approximately 10 minutes for the sample to come to bath temperature at temperatures up to 40°C and 15 minutes for temperatures up to 100°C.
7. Leaving the viscometer in the bath, apply suction to tube A and draw the liquid slightly above mark C (this can best be accomplished using a pipet bulb on tube A).
8. To measure the efflux time, allow the liquid sample to flow freely down past mark C, measuring the time for the meniscus to pass from mark C to mark E.
9. A check run may be made by repeating steps 7 and 8. Do a minimum of three measurements of efflux time.
10. Calculate the kinematic viscosity in mm²/s (cSt) of the sample by multiplying the efflux time in seconds by the viscometer constant.



7

© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory

- Each group takes the viscosity of two solutions at three temperatures.
- Data are shared with all lab classmates (Google Forms)
- See table for assignments (also given in the lab handout)

Spring 2016		
station	solutions to measure (wt% sugar)	
2	10.0	30.0
4	30.0	50.0
6	mystery	10.0
7	mystery	50.0
8	mystery	30.0
9	50.0	mystery
10	10.0	50.0

8

© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory

Lab: Measure Fluid Viscosity

- Measure densities, efflux times of sugar solutions at various concentrations, temperatures
- Calculate viscosities; enter into Google Form
- Determine experimental uncertainty
- Determine best-fit line at three temperatures; plot on one graph; compare to literature
- Identify concentration of mystery solution (with uncertainty)

9

© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory

Class: In your submissions to the Google Form:

- Do not write 30%; just write 30
- Do not write the units in your entries; it's going directly to Excel.
- If you want to use scientific notation, please write 3.14e3 for 3140 (do not write 3.14*10^3). It's going directly to Excel.

10

© Faith A. Morrison, Michigan Tech U.

CM3215 Fundamentals of Chemical Engineering Laboratory**Prelab assignment (due Monday):**

- Obtain the literature values of solution densities and viscosities before lab; this must be in your notebook for Monday's prelab check.
- **Please do the pre-lab**—this is to enhance your laboratory experience
- Note the references where you found your values.
- See instructions for complete prelab assignment.

11

© Faith A. Morrison, Michigan Tech U.