

Experiment 2: RTD Sensor

(tbc 1/14/2012)


Main Task:

Build an RTD circuit to measure temperature (to within $\pm 2^{\circ}\text{C}$).

Name: _____ Date _____

1. RTD circuit _____
2. Labview data acquisition program _____
3. Temperature Calibration Curve _____
4. Temperature test _____

1. The RTD Circuit.



- The 5 volts source can be obtained from the DIGITAL side of the USB-6008 device
- Use $R = 100 \Omega$ and connect only different colored wires of the RTD in the circuit. (Use a circuit breadboard to implement the circuit.)

2. The Labview program.

- Include **[Express]→[Execution Control]→[While loop]**. (Include everything below into the **[while loop]** window.)
- Import **[Express]→[Input]→[DAQ Assistant]** block. Choose **[analog input] → [voltage] → [ai0]**

- **Max** = **0.5 volts**
- **Min** = **0 volts**
- **acquisition mode** = **continuous**
- **samples to read** = **10**
- **rate** = **100 Hz**

- Import **[Express]→[Signal Analysis]→[Filter]** block.

- **Infinite Impulse Filter** = **selected**
- **Type** = **Lowpass**
- **Topology** = **Butterworth**
- **Order** = **3**
- **Cutoff Frequency (Hz)** = **1**

- d. Import **[Express]→[Signal Manipulation]→[from Dynamic Data]** block.
Choose the “single scalar”.
- e. Import **[Express]→[Signal Manipulation]→[to Dynamic Data]** block.
Choose “single scalar”.
- f. Import **[Express]→[Arithmetic & Comparison]→[Formula]**.
 - **Input X1, label** = v
 - **Formula** = v
- g. Include the following control blocks: **[Waveform Chart]**, **[Indicator]**.
- h. Match the wiring as shown in Figure 1.




Figure 1. Labview program.

3. Temperature Calibration Curve.

- a. Using a digital thermocouple, record the voltage readings corresponding to temperature readings at approximately 10°C intervals. (see Table 2 as example)

Table 1.

Temperature (°C)	Voltage (volts)
0	
25	
50	
75	
100	

- b. Obtain a curve fit (2nd order polynomial fit) of temperature as a function of voltage using Excel.

Formula:

- c. Modify the entry in the “**Formula**” block using the conversion formula obtained in step 2.
- d. Test the obtained **RTD VI**.