

STANDARD OPERATING PROCEDURES

***** ZEISS EM 10 *****

I. CONFIRM THAT THE FOLLOWING ARE TRUE:

- High voltage selection is at 80 kV
- Vacuum indicator is $< 10^{-5}$ Torr
- The green light is illuminated on the button labeled “HV ON”
- Set the magnification to read 100x.
- Toggle switch for the CCD camera is to the right (camera out)



II. SIGN INTO THE LOGBOOK

III. FILL THE LIQUID N₂ DEWAR FOR THE DECONTAMINATOR *(Always wear safety glasses when handling liquid nitrogen.) Be certain that the Shield Master is covering the viewing chamber glass.*




Check the status of beam before loading specimen (IV-V)

IV. TURN ON THE HIGH VOLTAGE AND HEAT THE FILAMENT:

- Depress the green “HV ON” button; the green light will go off and the red light will become illuminated on the button labeled “HV OFF”
- Wait 10 seconds for HT to stabilize.**
- Depress the button labeled “FIL”; the red light on this button will also become illuminated, and the meter that indicates the filament heating will jump up to about 1.0 V and settle to 0.6 V.
- Wait a few seconds to stabilize.
- Adjust the Condenser 2 knob to make a focused image of the illuminated spot
- If you cannot find the beam, Check:
 - CCD camera is out (switch to right),
 - Specimen not blocking beam.
 - Objective aperture out.

Note: If objective aperture is withdrawn beyond the open position, air leaks into the column. This is shown by a suddenly very weird beam image and very poor vacuum. This situation is to be avoided.

- With focused spot, increase magnification to 6,300x.
- Turn the filament heater knob clockwise to get an elliptical “eye” shape.**
- Use the “Y” and “X” Beam Alignment knobs to position the focused under-saturated illumination spot at a convenient position on the viewing screen;
- The focused spot should have an elliptical “eye” shape.
- If needed, adjust the C₂-STIGM. “Wheels” to get the sharpest image of filament.
- Decrease the magnification to lowest, 100x.
- Spread beam clockwise to fill viewing screen. Now you can load a specimen.



V. TURN OFF THE FILAMENT AND THE HIGH VOLTAGE


- Decrease (turn counterclockwise) the Filament Heating knob [ON THE LOWER LEFT control panel], all the way to the “stop”; the meter should read ~0.7.
- Press and release the button labeled “FIL”; the red light on this button will go off as the button pops out, and the meter that indicates the filament heating will drop to 0.0

- c. Press and release the button labeled “**HV OFF**”; the red light on this button will go off as the button pops out, and the green light will become illuminated on the “**HV ON**” button



VI. REMOVE SPECIMEN HOLDER AND INSERT YOUR SPECIMEN

- a. **Zero the specimen translators to 500, 500**
- b. Make sure filament and high voltage are **OFF**.
- c. Pull out specimen rod to forward position, turn counter-clockwise from 6:00 to 12:00 and pull out until it engages in stop position.
- d. Turn large knurled wheel fully counter-clockwise.
- e. Screw key into slide, remove slide by gently pulling on key, and place slide in slide support. **Do not touch slide or specimen holder with fingers.**
- f. Place cartridge upright and support with holder on slide support.
- g. Remove screw cap with tool.
- h. Place grid in cartridge, specimen side down.
- i. Replace screw cap. Watch bottom edge of cap move. Turn only until it stops. **Do not over-tighten or force cap.**
- j. Swing aside slide support and lower cartridge.
- k. Remove cap tool.
- l. Holding slide key, insert slide into air lock tube.
- m. Remove key.
- n. Turn large knurled wheel fully clockwise to seal air lock chamber.
- o. Press red button to release specimen rod and push in rod to stop position.
- p. Turn rod clockwise in milled-out ring **to 3:00 position. Rod will be in horizontal position. Red pre-pump light comes on. If light does not turn on or it red light does not turn off-- get Help - DO NOT PROCEED.**
- q. **When red light below meter for high vacuum/beam current goes out,** turn rod to 6:00 position.
- r. Press red button on specimen rod and push into intermediate stop. Vacuum level may jump a bit.
- s. Press red button again and gently push rod, 2 more times, to insert cartridge into stage.



VII. TURN ON THE HIGH VOLTAGE AND THE FILAMENT--Refer to the instructions above. **For viewing you sample heat the filament no more than a meter reading of 1.0 V; this is saturation for this filament.**

VIII. ENDING YOUR SESSION ON THE MICROSCOPE

- a. **Toggle switch for the CCD camera to the right (camera out)**
- b. Turn off the filament and the high voltage – see item V. above
- c. Remove your specimen from holder and put holder back into microscope.
- d. As courtesy--turn mag to 100x and spread C2 course fully clockwise to spread beam.
- e. Enter the filament time in the log book
- f. Write down any problems you found in the logbook and email problems to [manager](#).
- g. **Turn off the light switches on the right-hand control panel, and turn off the room lights as you leave the room.**

Contact Manager with any problems.

STANDARD OPERATING PROCEDURES

*****ZEISS EM 10*****

- **When leaving the room, look back at the TEM**
 - a. Is the camera out, to the right?**
 - b. Are the desk lights off?**
 - c. Is the room light off?**